

Fixed Depth Trepine

Nanoedge “Fixed Depth Trepine” are ready-to-use circular blades designed for accurate fixed depth cutting of corneal donor buttons or the recipient eye during corneal transplant surgery.



Product Description

Product Code	Description	Qty. per box (Sterile)
MIPL/C6/3.00	3.00 mm Diameter	1 Pc.
MIPL/C6/5.00	5.00 mm Diameter	1 Pc.
MIPL/C6/7.00	7.00 mm Diameter	1 Pc.
MIPL/C6/7.50	7.50 mm Diameter	1 Pc.
MIPL/C6/8.00	8.00 mm Diameter	1 Pc.
MIPL/C6/8.50	8.50 mm Diameter	1 Pc.
MIPL/C6/9.00	9.00 mm Diameter	1 Pc.

Each size is available in 150 µm, 200 µm, 250 µm, & 300 µm blade height.

Features

- An ultrasharp thin profile cutting edge minimizes damage to the endothelial cells
- Designed for a precise and symmetrical cut of the corneal tissue
- See through design helps in accurate positioning
- Can be used with or without an endothelial punch
- Innovative protective packaging to ensure safety of the cutting edge
- EO Sterile

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)