




Rose Touch (Rose Bengal Strips)

Rose Bengal is an ocular stain that has affinity for degenerating epithelium.



Product Description

Product Code	Description	Qty. per box
MIPL/A2	Rose Bengal Strips	100 strips (Sterile)

Indications

- It stains devitalized or degenerated epithelial cells as well as mucous filaments

Specifications

- Available as sterile strips in individual pouches
- Each strip is impregnated with 1.50 mg of Rose Bengal
- Innovative peel open packing with transparent film on one side

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)