



## Teflon Block

The corneo sclera tissue is placed over the Teflon block and the donor button then punched out.



## Product Description

Product Code	Description	Qty. per box (Sterile)
MIPL/D12	Teflon Block	10 pc. (Sterile)

## Features

- Aspheric curvature to match corneal curvature
- Central hole helps in accurate positioning
- Can be used with or without an endothelial punch
- EO Sterile
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)