


Capsule Polisher

Disposable ophthalmic cannula used for removing debris from the posterior, anterior and equatorial regions of capsular bag.

Micro etched tips are slightly roughened for gentle polishing.



Product Description

Product Code	Descripton with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E12/A1	Capsule Polisher (Kratz), 21G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A2	Capsule Polisher (Kratz) - 23G	22 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A4	Capsule Polisher (Kratz) - 25G	25 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A5	Capsule Polisher (Kratz) - 27G	25 mm	7 mm	10 Pcs. (Sterile)
MIPL/E12/A3	Capsule Polisher (Pearce) Olive Tip - 25G	25 mm	7 mm	10 Pcs. (Sterile)



Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)