

## CAM Vision Stimulator

This instrument is designed for the treatment of lazy/amblyopic eye based on the concept of active and controlled simulation.



### Product Description

Product Code	Description	Qty. per box
MIPL/M5	Cam Vision Stimulator	1 pc.

### Features

- Spatially oriented stimulus
- All-day passive occlusion replaced by seven-minute active stimulation
- Compact, portable model which can be given to the patient saving time and helping in regular exercise
- Size : 21 x 14 x 11 cm
- Weight: 1 Kg
- 6 Discs provided with the instrument
- CE Certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: [sales@madhuinstruments.com](mailto:sales@madhuinstruments.com) | [www.madhuinstruments.com](http://www.madhuinstruments.com)

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)