

Coaxial Cannula

Disposable ophthalmic cannula used for irrigation and aspiration of the anterior chamber.



Product Description

Product Code	Descripton with Gauge	Gauge	Port Size	Qty. per box
MIPL/E19/A2/01	Coaxial I&A Handpiece (Straight) - Disposable	21G	0.40mm	1 Pc. (Sterile)
MIPL/E19/A2/02	Coaxial I&A Handpiece (45° Angle) - Disposable	21G	0.40mm	1 Pc. (Sterile)
MIPL/E19/A3/01	Coaxial I&A Handpiece (Straight) - Disposable	21G	0.30mm	1 Pc. (Sterile)
MIPL/E19/A3/02	Coaxial I&A Handpiece (45° Angle) - Disposable	21G	0.30mm	1 Pc. (Sterile)



Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)