



Flute Handle (Charles)

This is a simple instrument for passive aspiration of intraocular fluids.



Product Description

Product Code	Description	Qty. per box
MIPL/N8	Flute Handle (Charles)	1 pc.

Features

- Supplied unsterile
- User friendly design/ergonomic grip
- CE Certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)