



Peroperative Keratoscope

Peroperative keratoscope is used during an ophthalmic surgical procedure to judge the astigmatism induced in the cornea while putting sutures.



Product Description

Product Code	Description	Qty. per pack
MIPL/N5	Peroperative Keratoscope	1 Pc. (Unsterile)

Features

- This is a mechanical device made out of anodized aluminum tube.
- The inner wall of the tube has alternate black and white rings of different thickness.
- The reflection of these rings on the patient's cornea helps in analyzing the astigmatism induced due to tight and loose sutures.
- Can be resterilised by autoclave or ETO
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Cleaning and Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)