

Phaco Practice Eye Fixation Head

This artificial head is used for fixation of the phaco practice eyes and give a more natural feel of the human face to the user.



Product Description

Product Code	Description	Qty. per box
MIPL/M6	Phaco Practice Eye Fixation Head	1 Pc.

Features

- Simulates facial contours
- Includes a fluid collection tray and a drainage tube
- Easy to clean
- Light weight – easy to carry
- Very economically priced
- CE Certified



Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)