

Three Mirror Goniolens

This gonio lens is used for viewing the anterior chamber angle and central & peripheral retina.



Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/I3	Three Mirror Goniolens	59°/66°/73°	1 pc.

Features

- This lens has three mirrors angled at 59°, 66° and 73° and arranged at 120° intervals to permit viewing of the fundus and anterior chamber angle.
- The posterior pole is also observable through the central axis of the lens.
- Made in glass.
- Contact medium required.
- CE Certified.

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)