

# Bimanual Cannula

Disposable ophthalmic cannula used for irrigation and aspiration of the anterior chamber.

## Product Description



Bimanual I & A Cannula - Disposable

Product Code	Description	Gauge	Port Size	Qty. per box
<b>MIPL/E15/DP/21WP</b>	<b>Bimanual Cannula Handpiece Set with Polisher – Disposable</b>			<b>1 Set (Sterile)</b>
MIPL/E15/A8/01	Bimanual Irrigation Handpiece – Disposable	21G	0.50mm	1 Pc. (Sterile)
MIPL/E15/A8/02	Bimanual Aspiration Handpiece with polisher – Disposable	21G	0.40mm	1 Pc. (Sterile)
<b>MIPL/E15/DP/22WP</b>	<b>Bimanual Cannula Handpiece Set with Polisher – Disposable</b>			<b>1 Set (Sterile)</b>
MIPL/E15/A9/01	Bimanual Irrigation Handpiece – Disposable	22G	0.45mm	1 Pc. (Sterile)
MIPL/E15/A9/02	Bimanual Aspiration Handpiece with polisher – Disposable	22G	0.35mm	1 Pc. (Sterile)
<b>MIPL/E15/DP/23WP</b>	<b>Bimanual Cannula Handpiece Set with Polisher – Disposable</b>			<b>1 Set (Sterile)</b>
MIPL/E15/A10/01	Bimanual Irrigation Handpiece – Disposable	23G	0.40mm	1 Pc. (Sterile)
MIPL/E15/A10/02	Bimanual Aspiration Handpiece with polisher – Disposable	23G	0.30mm	1 Pc. (Sterile)



Bimanual I & A Cannula - Titanium

Product Code	Description	Gauge	Port Size	Qty. per pack
MIPL/E15/A6/01	Bimanual Irrigation Handpiece - Ti	21G	0.50mm	1 Pc.
MIPL/E15/A6/02	Bimanual Aspiration Handpiece with polisher - Ti	21G	0.40mm	1 Pc.
MIPL/E15/A7/01	Bimanual Irrigation Handpiece - Ti	22G	0.45mm	1 Pc.
MIPL/E15/A7/02	Bimanual Aspiration Handpiece with polisher - Ti	22G	0.35mm	1 Pc.
MIPL/E15/A5/01	Bimanual Irrigation Handpiece - Ti	23G	0.40mm	1 Pc.
MIPL/E15/A5/02	Bimanual Aspiration Handpiece with polisher - Ti	23G	0.30mm	1 Pc.



## Features

- Ergonomic grip
- Thin wall cannula used for maximum flow
- Standard male/female luer connector

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)