



Tear Touch (Schirmer Tear Test Strips)

This diagnostic strip made from the Whatman No 41 filter paper measures the amount of tear production.



Product Description

Product Code	Description	Qty. per box
MIPL/A1	Schirmer Tear Test Strips	100 strips (Sterile)

Indications

- Used for measurement of Tear Production

Specifications

- Each strip has a printed mm scale so that they can be easily read and recorded
- Available as sterile strips in individual pouches
- Innovative peel open packing with transparent film on one side
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)