

## 20D Aspheric Lens

This slit lamp lens is used for general examination of the fundus.



### Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/K1	20D Aspheric Lens	50°	1 pc.

### Features

- This Slit lamp lens offers an ideal combination magnification and field of view.
- It also provides an extended working distance from cornea.
- Superior glass optics offers high clarity images.
- Made in glass.
- Aspheric Optics

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

### Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: [sales@madhuinstruments.com](mailto:sales@madhuinstruments.com) | [www.madhuinstruments.com](http://www.madhuinstruments.com)

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)