



## Endothelial Punch

This mechanical device gives an excellent donor corneal button punched out from the endothelial side of corneoscleral button.



### Product Description

Product Code	Description	Qty. per box
MIPL/N1	Endothelial Punch	1 Pc.

### Features

- Universal - Any size of trephine from 2 mm to 11 mm can be fitted
- Precentered - The trephine and the teflon block are automatically centered for all sizes of trephines
- Spring Recoil - After completion of the cut, the trephine retracts back with an inbuilt spring system
- Sterilization - Endothelial Punch can be sterilized by autoclaving or EO gas
- Each set consists of two plungers

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: [sales@madhuinstruments.com](mailto:sales@madhuinstruments.com) | [www.madhuinstruments.com](http://www.madhuinstruments.com)

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)