




Irrigating Vitrectomy Lens Set

This set comprises of 3 lenses which have a provision to attach the infusion line for irrigation during a VR surgery.



Product Description

Product Code	Description	Field of View	Magnification	Contact Diameter	Qty. per box
MIPL/L2/BC	-90D/Biconcave		0.80 X (air filled)	10 mm	1 pc.
MIPL/L2/WF	Widefield		0.50 X (fluid filled)	10 mm	1 pc.
MIPL/L2/CF	Central Fundus		1.00 X	10 mm	1 pc.

Features

- High quality glass optics for clear images
- Lens description engraved on the handle for easy identification
- Flexibility to order individual lenses
- CE Certified
- Sterilization – can be Sterilized by Autoclaving or EO Gas

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)