

## YAG Iridectomy Lens

This lens is used for YAG laser iridectomy procedures.



### Product Description

Product Code	Description	Qty. per box
MIPL/J2	YAG Iridectomy Lens	1 Pc.

### Features

- This lens consists of a modified Goldmann type fundus lens with a 7.50 mm diameter +70 D plano-convex lens bonded to the anterior surface.
- The +70 D button lens helps in condensing the laser energy so that the procedure can be carried out at lower power settings.
- The anodized light metal rim gives a firm grip on the lens.
- Made in glass.

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

### Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: [sales@madhuinstruments.com](mailto:sales@madhuinstruments.com) | [www.madhuinstruments.com](http://www.madhuinstruments.com)

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)